

MT500

MOBILE TABLE

The DEL MEDICAL MT500 mobile radiographic table is ideally suited for U-Arm systems located in imaging centers, orthopedic clinics and private practices. It can also be paired with an overhead tube crane or floor mounted tube stand to offer an economical table option. It features a heavy-duty steel frame for maximum stability and a stain-resistant fiber-resin tabletop. The large-diameter casters provide smooth movement, swivel 360 degrees and come standard with easy, yet reliable, foot locks.

Key Highlights

4 Durable Swivel Casters with Foot Locks

All Steel Structure Welded Base

- Provides maximum strength and reliability

500 lbs. Patient Weight Capacity

Smooth, Stain Resistant Fiber Resin Top






Technical Specifications

* All specifications subject to change without notice

MT500 Mobile Table

Load capacity:	500 lbs. (227 kg)
Tabletop length:	27.9" (70.9 cm)
Tabletop width:	26" (66.1 cm)
Tabletop height:	27.9" (70.9 cm)
Surface to cassette tray film plane:	3.2" (8.0 cm)
Weight without cabinet:	164 lbs. (74 kg)
Weight with cabinet:	201 lbs. (91.2 kg)
Options:	Grid cabinet with tray (focal spot travel length of 43.0" (109.22 cm)) / Table-side patient handgrips / Lateral cassette holder / Rotating tray / Compression band



DEL MEDICAL



#9000-MT500-DS 09.20 Rev. B

New York

28 Calvert Street
Harrison, NY 10528
(Tel) 800.261.9808
(Fax) 914.835.6111

Illinois

241 Covington Drive
Bloomington, IL 60108
(Tel) 800.800.6006
(Fax) 847.288.7011

www.delmedical.com