DEL MEDICAL



VS100

VERTICAL WALL STAND

A counterbalanced wall stand with vertical electric locks that allow for easy positioning. The extensive vertical travel allows for a complete range of exams from skull to weight bearing knee and ankle exams.

Key Highlights


Electromagnetic Locks

- Conveniently placed for easy vertical positioning

Dual Counterweight Cables for Added Safety

Vertical Travel Range of 58.4" (148.3 cm)

Compatible with ISO 4090 (cassette-sized) detectors as well as standard CR and analog film cassettes of various sizes




VS100 Wall Stand	
Column height:	85.1"(216.1 cm)
Receptor center height max:	72.4" (183.9cm)
Receptor center height min:	14"(35.6 cm)
Receptor travel:	58.4" (148.3 cm)
Wall to receptor distance:	10.85" (27.6 cm)
Front panel aluminum equivalency:	0.4 mm
Weight:	280 lbs. (127 kg)
Power requirements:	24 V DC at 0.5 amps


Certifications

ETL listed

VS100 Options:

- 14" x 17" DR tray
- 17" x 17" DR tray
- Standard tray for smaller sized CR or analog cassette
- Size-sensing cassette tray with PBL interface
- Ion Chamber (AEC)
- Variety of grid options







(Fax) 847.288.7011