

FMT

FLOOR MOUNTED TUBE STAND

The FMT tube stand is easy to use and features a continuous digital angulation display. SID is visible when the system is positioned at 40 or 72 in (102 or 183 cm) to the wallstand or 40 in (102 cm) to the table. Conveniently positioned finger-tip control buttons, and powerful magnetic locks ensure steady, accurate tube head positioning. In addition, it only requires a ceiling height of 7 ft 3 in (221 cm) making it one of the lowest ceiling height requirements in the industry.

The FMT System was developed to keep costs low while also producing the highest quality diagnostic images. Built in the USA, it will deliver many years of continuous use and exceptional user-satisfaction.



Key Highlights

Column Rotation of +/-180°

Manual Collimator (Automatic Available)

Electromagnetic Brake System

- For vertical and longitudinal travel

Standard 10 ft Floor Plate

- Optional 8 or 12 ft floor plates available

11 in Transverse Travel

- Access to full table top coverage

Minimum Focal Spot Positioning at 10 in


- Perfect for weight bearing exams

Technical Specifications

* All specifications subject to change without notice

FMT Tube Stand

Minimum Ceiling Height:	87 in (221 cm)
Column Height:	85.5 in (217 cm)
Vertical Travel:	64.5 in (164 cm) total
Longitudinal Travel:	96 in (244 cm)
Transverse Travel:	11 in (28 cm)
Floor-to-Focal Spot (with Standard Platform Mount):	Minimum 10 in (25.4 cm), Maximum 74.5" (190 cm)
Rail Longitudinal Length:	10 ft (3 m) / Optional 8 ft (2.4 m) or 12 ft (3.7 m)
Column Rotation:	180° where installation permits
Tube Rotation:	±180° detents for lateral decubitus projections
Base Assembly to Rear Wall:	8 in (20.3 cm)
Tube, Collimator, and Cable Support Capacity:	88 lbs (40 kg)
Weight (Excluding X-ray Tube, Collimator, and Counterweights):	200 lbs (91 kg)
Power Requirements:	120 / 230-240 VAC 50 / 60 Hz



DEL MEDICAL



#9000-FMT-DS 02.21 Rev. C

New York

28 Calvert Street
 Harrison, NY 10528
 (Tel) 800.261.9808
 (Fax) 914.835.6111

Illinois

241 Covington Drive
 Bloomingdale, IL 60108
 (Tel) 800.800.6006
 (Fax) 847.288.7011

www.delmedical.com