


VT300 *Family*

Tilting Wall Stand w/Optional Vertical Tracking



Key Highlights

Receptor Tilt Range: +90° / -20°

Available with Cabinet Options

- Flex Tray, standard cassette tray, rotating tray, or 17 x 17 in (43 x 43 cm) fixed detector mounted inside

Electromagnetic Locks

- Conveniently placed control buttons are positioned on the side of the receptor cabinet allowing for easy maneuverability and vertical travel

Optional AutoTracking (FMT18T and OTC18T/S)

Vertical Travel Range of 56 in (142 cm)

- Provides patient comfort for standing or seated chest examinations

Flexible Installation

- Control buttons and load side can be changed to either side of the image receptor

Smooth, Counterbalanced Movements

- Accommodating to all users regardless of size and strength

Motorized Vertical Tracking (VT300T)


VT300 Family

TECHNICAL SPECIFICATIONS

Receptor Tilt Range:	+90°/-20°
Vertical Travel Range (Vertical Receptor):	15.5 in (39.4 cm) to 71.1 in (180.6 cm)
Vertical Travel Range (+ 90° Receptor):	25.5 in (64.8 cm)
Height:	VT300: 82.24 in (208.9 cm) VT300T: 86.5 in (218.9 cm)
Depth:	30.17 in (76.6 cm) vertical / 40.1 in (101.9 cm) horizontal
Receptor Width:	26.5" (67.2 cm)
Weight:	765.5 lbs. (347.2 kg) total 368.5 lbs. (167 kg) without counterweights
Power Requirements:	24 VDC 4 A (from DM-PSI)

Environmental Requirements

	Operating	Transit/Storage
Temperature:	50°F - 95°F (10°C - 35°C)	-40°F to 158°F (-40°C to 70°C)
Relative Humidity Limits:	10% to 80% non-condensing	10% to 100%
Atmospheric Limits:	14.5 inHg to 30.74 inHg / 800 hPa to 1060 hPa	



DEL MEDICAL



#9000-VT300 Rev. 3

New York

28 Calvert Street
Harrison, NY 10528
(Tel) 800.261.9808
(Fax) 914.835.6111

Illinois

241 Covington Drive
Bloomington, IL 60108
(Tel) 800.800.6006
(Fax) 847.288.7011

www.delmedical.com